

## **2. Solid-Liquid separation method- Filtration**

**Aim** : Is to determine the flow rate of the tissue homogenate solution through filter media.

### **Principle:**

- Filtration is ubiquitous in solid liquid separation of suspension mixer.
- Liquid to be clarified is passed through a porous barrier.
- Large particles are retained at the filter surface or within the depth of filter medium; solvent and some small particles pass through the filter.
- Process is pressure driven and unlike sedimentation works well even when there is no density difference between the particulated and the suspending medium.
- When membranes tested using pure water most membranes give very high fluxes.
- But during actual filtration the flux is likely to be, very much lower.
- At the start of operation the flux is high and in some cases may approach the pure water flux.
- However, a rapid decline occurs over the first 20-30 min of operation followed by a gradient, but persistent, fall.
- This reduction in flux is caused by a combination of concentration polarization and fouling.

## Equipments and reagents:

1. Sheet of membrane to test
2. Ground cell Suspension which can be loaded with flat sheet of membrane.
3. Ancillary Equipment for operation of membrane module.
4. Timer.



## Determination of flow rate of homogenate solution through filter media:

- 1) Close all the inlet and outlets of permeate reservoir and Sample reservoirs of the Ancillary Filtration Unit
- 2) Fill tissue homogenate solution to the Sample reservoir.
- 3) Cut out a disc of membrane to fit into a dead-end filtration unit.
- 4) Wash thoroughly to remove preservatives and wetting agents and fix the disc to the Filter Membrane holder.
- 5) Load the tissue homogenate solution through the membrane by opening the inlet and outlet valves of the sample reservoir and switching on the pump.
- 6) Apply Fixed Pressure and measure the permeate flow rate by measuring volume filtered over period of 5, 10, 15, 20, 25 and 30 sec.
- 7) Plot a graph by plotting flow rate against time.

### **Results:**

- |  |             |
|--|-------------|
| 1. Flow rate of the distilled water after 2000 sec:  | 0.2 cc/sec  |
| 2. Flow rate of the tissue homogenate after 2000 sec | 0.07cc/sec  |
| 3. Difference in the flow rate:                      | 0.13 cc/sec |



OBSERVATIONS

Time (s)	Volume of sample collected [cc]	Flow rate [cc/s]
300	40	0.133
600	70	0.116
900	95	0.105
1200	110	0.092
1500	120	0.08
1800	130	0.072

Note: plot a graph of flow rate v/s time in sec both for membrane fouling and cake formation



